



Make the invisible visible



# Using the Discovery-Q to model liver fibrosis; quantitatively measure the effect of co-cultures on liver cell health and behavior

## Abstract

The Discovery-Q is a biosensor that measures the changes in mass and viscoelasticity of adherent cells in real-time over multiple days. These measurements relate to co-culture homeostasis and cell-to-cell tension cell-to-extracellular matrix tension. The Discovery-Q allows for hepatocyte and stellate co-cultures ratios mimicking physiological disease states. This modeling allows for generation of signatures of disease states and testing of the response of drugs on this model.


**“Based on data presented, we believe that this method would allow for detection of cellular changes following toxicant treatment; thus, becoming a real-life platform for drug detection and testing.”**

Rafal Witek, Director Advanced Cell Systems, Thermo Fisher Scientific.




Fibrosis co-cultures with hepatocyte and stellate ratios to simulate liver fibrosis. The Discovery-Q can quantify the effect of adding stellates to the culture and the resulting rigidity and tension in the cell layer mimicking liver disease states.

## Modeling co-culture and cellular interplay in disease progress



Hepatocyte & Stellate interaction



Monitoring Cell to cell effects



Monitor Real-time Cell Health



## Product information


Platform: Discovery-Q

- Compatible with adherent and semi-adherent cell lines
- Simultaneous real time readings of frequency and resistance (dissipation)
- Minimum cell number per well: 1,000
- Label free detection, **specialized media or serum not required**
- Network based system
- Remote operations on web browser
- Data is date and time stamped and downloadable as CSV files


Cell culture types tested: Primary cell lines (cryopreserved plateable hepatocytes, stellates and co-cultures), isolated cells from liver, breast, lung, umbilical, macrophage, and primary tumors

Cell culture lines tested (not a comprehensive list): HepG2, Sk-Hep-1, DH82, SK-Br-3, MDA-MB-231, HUVEC, BAE, BAEC, HL-60, HT-29, HMEC, NHBE, HMVEC-L, HMVEC-BL, HepRG, Hs578t, FaDu, MCF-7, MCF-12A,

For information regarding pricing of units or any other matters please contact: [info@invitrometrix.com](mailto:info@invitrometrix.com)



The Discovery-Q and consumable well plate



Invitro-Q units working in a cell culture incubator.



**Invitrometrix Corporation**  
116 John St. Ste. 340  
Lowell, MA 01854  
[www.invitrometrix.com](http://www.invitrometrix.com)

For Research Only:  
Not for use in diagnostic procedures

This information is subject to change without notice.  
© Invitrometrix Corporation 2019  
Published in the USA, July 2019  
AB-004